

Standard Specifications


| Model | Total Organic Carbon Analyzer TOC-300V | | | | | | | | | | | | | | |
|--|---|-----------------------|---------------|-------------------------------|-----------------------|---------------|-----|-----|----|----|-----|----|-----|-----|----|
| Samples | Water Samples (Drinking, Industrial Waste, Environmental) | | | | | | | | | | | | | | |
| Method | Oxidative Combustion/NDIR Detection | | | | | | | | | | | | | | |
| Furnace | max. 900℃ | | | | | | | | | | | | | | |
| Measurement Items | TC, IC, TOC, NPOC (optional TN) | | | | | | | | | | | | | | |
| Measurement Range | 0.1 to 1000mg/L (TC, IC). Auto-dilution for >300mg/L. | | | | | | | | | | | | | | |
| Sample Amount | 50 to 500μl | | | | | | | | | | | | | | |
| Measurement Time | <4min. (depends on sample amount and parameters) | | | | | | | | | | | | | | |
| Autosampler | 60 positions | | | | | | | | | | | | | | |
| Measurement Accuracy | <table><tr><th>Concentration TC·IC (mg/L)</th><th>Sample amount (μl)</th><th>RSD(%) n=5</th></tr><tr><td>1.0</td><td>100</td><td><8</td></tr><tr><td>10</td><td>100</td><td><5</td></tr><tr><td>100</td><td>100</td><td><3</td></tr></table> | | | Concentration TC·IC (mg/L) | Sample amount (μl) | RSD(%) n=5 | 1.0 | 100 | <8 | 10 | 100 | <5 | 100 | 100 | <3 |
| | Concentration TC·IC (mg/L) | Sample amount (μl) | RSD(%) n=5 | | | | | | | | | | | | |
| | 1.0 | 100 | <8 | | | | | | | | | | | | |
| | 10 | 100 | <5 | | | | | | | | | | | | |
| | 100 | 100 | <3 | | | | | | | | | | | | |
| TC: Potassium Hydrogen Phtalate Standard Solution IC: Na2CO3/NaHCO3 Standard Solution | | | | | | | | | | | | | | | |
| Gas | High purity air or standard air. Oxygen necessary for optional TN measurement (O2 >99.7%). | | | | | | | | | | | | | | |
| Power | Main Unit: AC100 to 240V, 500VA Autosampler: AC100 to 240V, 80VA | | | | | | | | | | | | | | |
| Dimensions | 750(W) x 530(D) x 1003(H)mm, 66kg | | | | | | | | | | | | | | |

Optional Unit

Nitrogen Detector ND-210

By attaching our well-proven ND-210 (developed for trace nitrogen analysis) to the TOC-300V system, total nitrogen and TOC can be measured simultaneously.

| | |
|------------|--|
| Method | Oxidative Combustion/Chemiluminescence Detection (Reduced Pressure Method) |
| Power | AC 100 to 240V |
| Dimensions | 220(W) x 375(D) x 500(H)mm |
| Weight | 22kg |



*The contents of this brochure are subject to change without notice.
*Company and product names contained herein are the trademarks or registered trademarks of the company concerned.

MITSUBISHI CHEMICAL ANALYTECH CO., LTD.

Instruments Division
370 Enzo, Chigasaki, Kanagawa 253-0084, Japan.
Tel: +81-467-86-3864 Fax: +81-467-86-3862
URL: <http://www.mccat.co.jp>

DIRECT INJECTION therefore TOUGH

TOC-300V

Total Organic Carbon Analyzer



Are your current TOC analyzers giving you the TRUE results?
Especially when it comes to samples with high amount of particulates?


We have developed a new TOC analyzer with DIRECT INJECTION METHOD.
This makes our instrument predominantly TOUGH against samples with particu-
lates. Our proficient combustion technology and stable detection technology will
also contribute to bring you the TRUE results.

Official Test Methods: ISO 8245 (EN1484), EPA 415.1, EPA 9060A,
Standard Methods 5310, ASTM D7573

Characteristics

TOUGH therefore FAST

Using the standard equipped autosampler and newly devel-
oped DIRECT INJECTION method, samples containing particu-
lates can be measured without any pre-treatment.
(patent pending)



TOUGH therefore ACCURATE

The system easily passes the so-called 'Cellulose Test' (ISO 8245, Annex B3). Stable recovery of organic carbon in sus-
pended solids (SS) can be achieved without any special treat-
ment.

1. Test Suspension
225mg/L Cellulose powder
(=100mg/L as TC)

2. Test Criteria
Recovery within ±10%
RSD < 10%

| Reps | Recovery |
|---------|----------|
| 1 | 95.8 |
| 2 | 96.7 |
| 3 | 95.3 |
| 4 | 96.5 |
| 5 | 96.7 |
| Average | 96.2 |
| RSD (%) | 0.64 |


TOUGH therefore SIMPLE

By simply opening the front door
of the main unit, reagents, com-
bustion tube, waste bottle, and
tubings are all easily accessible
for checking and maintenance.




TOUGH therefore CLEAN

The complete sample channel is VALVE FREE and the sample
never enters the syringe. You are free from worries about clog-
ging and cross contamination.



TOUGH therefore FRIENDLY

Simple configuration and icons of the software makes the
system easy to use even at first sight.
Our environmentally friendly design philosophy has resulted in
a small furnace with low electricity consumption.




TOUGH therefore CONVENIENT

Just set your sample and start!
The auto-dilution function will
automatically dilute your
unknown sample in case the
initial result exceeds the cali-
bration range.



TOUGH and FAST for NPOC Measurements

The TOC-300V has achieved the FASTEST NPOC measurement in the industry
with the auto acid addition and 3-vial simultaneous pre-treatment.



Applications

| Sample | Mode | Reps | | | Average (mg/L) | RSD (%) |
|------------------------------|------|-------|-------|-------|-------------------|------------|
| | | 1 | 2 | 3 | | |
| River water A | TOC | 3.256 | 3.354 | 3.211 | 3.273 | 2.23 |
| River water B | TOC | 1.266 | 1.232 | 1.388 | 1.295 | 4.99 |
| River water C | NPOC | 2.187 | 2.283 | 2.214 | 2.228 | 1.93 |
| River water D | NPOC | 3.541 | 3.721 | 3.748 | 3.670 | 2.80 |
| Industrial waste water A | NPOC | 4.357 | 4.620 | 4.452 | 4.476 | 2.13 |
| Industrial waste water B | NPOC | 1.007 | 1.115 | 0.990 | 1.038 | 4.92 |
| Potassium Hydrogen Phthalate | TOC | 98.54 | 99.12 | 99.31 | 98.99 | 0.41 |

Measurement Principles

■ Oxidative Combustion-NDIR Detection Method

Organic compounds in the sample are combusted at high tem-
perature in high purity air in the presence of oxidation catalyst.
The carbon dioxide resulting from the combustion is detected by
NDIR to measure TOC.
TOC = TC - IC

■ What is NDIR (Non-dispersive Infrared Sensor)?

CO₂ molecules selectively absorb infrared rays with 4.26
micron wavelength. This absorption is proportional to the con-
centration of CO₂. This characteristic is used to quantify CO₂
with NDIR.

TC (Total Carbon) Measurement

All carbon (organic and inorganic) in the sample is either combusted or decomposed to carbon dioxide (CO₂) by injecting the
sample into a high temperature combustion tube filled with oxidation catalyst. The CO₂ is detected by NDIR and detection
signal is converted to a peak. The TC concentration in the sample can be obtained by correlating the peak area to a calibration
curve created from standard solutions.

IC (Inorganic Carbon) Measurement

Sample is injected into an IC reaction chamber which is filled with IC reaction solution (25% Phosphoric Acid). Due to the acidic
environment, only the inorganic carbon in the sample is converted to CO₂. The CO₂ is again detected by NDIR and IC concentra-
tion is calculated in the same way as TC concentration.

NPOC (Non-Purgeable Organic Carbon) Measurement

Obtaining TOC by subtracting IC from TC can result in high error, especially for samples containing high amount of IC and
relatively low amount of TOC. To avoid this error, NPOC is often used as a substitute and in many cases regarded as
'equivalent' to TOC. To run an NPOC measurement, sample is first treated with acid, then sparged to remove IC, followed by
a TC measurement.